




User Manual



VENUS

Walker

Product Code: WA03013

CareCo, Hubert Road, Brentwood, Essex, CM14 4JE

General Warning

Thank you for purchasing the Tuni Venus walker.

1. DO NOT install or use this walker without first reading and understanding these instructions. If you are unable to understand the warnings, cautions or instructions in this manual, please contact a healthcare professional, dealer or technical personnel before attempting to install this walker – otherwise, injury or damage may occur.
2. Always make sure the handle height on both sides is the same and the knob there is tightened. Otherwise, it may cause injury.
3. Please do not open this walker with great force as this may cause damage to parts.
4. Do not lean too far over the front of the walker to prevent instability and the risk of fall.
5. Do not use this walker if the weight of user exceeds the maximum. Weight capacity of this walker is 136 kg.
6. Always check the condition of parts to ensure safety.

Features for this product

1. The Tuni Venus is modern looking, user-friendly, improves stability, strong and durable.
2. It folds compact for travel and easy storage.
3. The walker is height adjustable.
4. Front wheels are attached for easier walking.

Specification

Folding size:

Width: 85.5 cm


Depth: 11.5 cm

Height: 85-95 cm

Height adjustable from 85 cm to 95 cm




Weight capacity: 136 kg (21st 5lbs)

Instructions



Assembly

- A. Install side frame into foot tubing.
- B. The hole on side frame is matched with that on foot tubing (handle height on both sides is the same) and then tighten the knob.
- C. Grasp the handle and open the walker to right position an audible "click" is heard as shown in figure below.
- D. Assemble bag tubing to each side of foot tubing and then hang the shopping bag.



Handle height adjustment

Loosen the knob on both sides and adjust the side frame to desired height, tighten the knob on condition that the handle height on both sides is the same.

Walking

Grasp handgrip upward with gentle force, push forward and then put down the walker. Support the body with handgrip and walk forward. Repeat this action to walk forwards, and keeping balance during usage.

Folding

Grasp the handgrip and pull the lever upward as shown in fig. 1 (Do not use lever during usage), and fold the handgrip inward as fig. 2 shows.

Maintenance

1. If something goes wrong with the walker, do not disassemble as this may cause damage to other parts. Please contact dealer or manufacturer.
2. Keep the walker away from acidic substances. Clean the surface with a mild soap detergent.